



TP IV LABORATORY REPORT

---

# Eda napari: implementation of a napari plugin

a plugin for visualising framerates and time from timelapse images

---

Steven Brown  
steven.brown@epfl.ch  
Supervisor: Willi Stepp

Laboratory of Experimental Biophysics, EPFL  
June 15, 2022

## Contents

<b>1</b>	<b>Introduction</b>	<b>2</b>
<b>2</b>	<b>Napari</b>	<b>2</b>
<b>3</b>	<b>Creating a plugin</b>	<b>3</b>
<b>4</b>	<b>Interactions between plugin and napari</b>	<b>4</b>
<b>5</b>	<b>Eda-napari</b>	<b>5</b>
5.1	Frame rate widget . . . . .	5
5.2	Time scroller widget . . . . .	6
5.3	Results . . . . .	7
5.4	Limitations and Improvements . . . . .	9
5.5	Automatic testing . . . . .	9
5.6	Issues . . . . .	10
5.7	Installing the eda napari plugin . . . . .	10
<b>6</b>	<b>Conclusion</b>	<b>10</b>
<b>7</b>	<b>Annexes</b>	<b>12</b>
7.1	Project structure . . . . .	12
7.2	The main widget python code . . . . .	12
7.3	Debugging . . . . .	20

## 1 Introduction

During biological observations numerous images are captured and the data is stored. The user must then choose a programme for viewing and analysing the data. Napari, a python based open source image viewer, is adapted for biology and facilitates the addition of custom plugins [1]. The aim of this project is to programme, following the napari architecture, a plugin for visualising metadata from image sequences.

More specifically, a plugin that is designed for images captured with an event driven acquisition algorithm. Capturing dynamical processes in biology can be carried out with the aid of time-lapse videos. However, these processes can take a certain amount of time to occur and the dynamics of these processes that take place may not be linear in time. This motivates the usage of event driven image acquisition. When activity is detected, then images are taken with a higher frame rate. This has several advantages. Firstly, large data is stored only for relevant time frames. If fluorescent imaging methods are used, this will also have the effect of minimising fluorophore bleaching, and thus enabling longer observations times [2].

The **eda-napari** plugin will be designed for visualising variable frame rates from image sequences. The programme has direct applications at the EPFL experimental laboratory of biophysics. In the laboratory, a neural network for detecting activity and adapting frame rates is used with an iSIM microscope for observing mitochondrial fission [2]. The plugin will allow visualising such data, acquired with adaptive capture times.

## 2 Napari

**Napari** is a multi-dimensional image viewing and analysis software [3]. It is an open and free to use software, that started back in 2018. Napari filled the need of a python based n-dimensional viewing software. Previously, users needed time consuming procedures to navigate between non-python based analysis software such as ImageJ and useful python packages [1]. Python has a great deal scientific packages that can be applied to image analysis, for instance deep learning packages scikit-learn, TensorFlow and PyTorch. Napari, being python based, can therefore easily integrate machine learning analysis tools. Napari's graphical interface is written with Qt. All sort of user interaction can thus be implemented within the Qt framework.

Napari itself can directly be used as an image viewer. Images, for example tiff files, with several layers can be displayed and napari holds the basic tools for colouring the image layers. Napari's main display consists of a main window, a scroll bar and side docks (Figure 1). The programme is especially created for allowing the user to install additional plugins. These plugins are simple to install and add image analysis tools, for example image segmentation, tracking or specific file readers.

For sharing plugins, napari is linked to **napari-hub**. Napari-hub allows users to have access to a whole environment of plugins. The plugins can be installed into napari and with the possibility of giving feedback. Napari-hub facilitates the search for plugins and improves communication between plugin creators and users.




Figure 1: Example of napari’s viewer display, visualising 3D + 2 channels cell images [4]

### 3 Creating a plugin

As previously said, it is possible to install plugins to napari and to add extra features to the base image viewer. This section will describe the main points for creating a plugin and is based on napari’s own tutorial<sup>1</sup>. For napari to detect and include a plugin, a specific minimum structure needs to be followed (Figure 2). A napari plugin can be seen simply as a python package and must have certain required files for installation.




Figure 2: Starting project directory tree for creating a npe2 napari plugin [4]

The *pyproject.toml* and *setup.cfg* files are specific to the package installation. Pyproject.toml uses setuptools to package the plugin, whereas setup.cfg defines package metadata, install requirements and sets a napari manifest as entry point to allow napari to find the plugin (Code 1). Napari has a second generation plugin engine npe2, with the old one soon to be deprecated<sup>2</sup>. With npe2, a *manifest (.yaml)* file must be defined. This file allows napari to discover the plugin and sets the different contributions. The contributions can be widgets, readers or simply commands. For clarity, plugin classes are often separated into multiple .py files for different types of contributions. For our example, we will consider a *\_widget.py* file (Code 1). The widget file would contain classes that define widgets, for example a simple button with interactions.

<sup>1</sup>Your First Plugin: [https://napari.org/plugins/first\\_plugin.html](https://napari.org/plugins/first_plugin.html)

<sup>2</sup>napari-plugin-engine was the first generation version. It used hook implementations "@napari\_hook\_implementation" to allow napari to detect the plugin.

Code 1: Example code of pyproject.toml, napari.yaml and setup.cfg, inspired from napari [4]  
**/pyproject.toml**

```
1 [build-system]
2 requires = ["setuptools", "wheel"]
3 build-backend = "setuptools.build_meta"
```

**/napari.yaml (manifest)**

```
1 name: plugin-name
2 contributions:
3   commands:
4     - id: plugin-name.command_name
5       title: command title
6       python_name: plugin_name._widget:widget_class_name #path of class or
7         function
8   widgets:
9     - command: plugin-name.command_name
10       display_name: Title in napari tab window
11       autogenerate: true
```

**/setup.cfg**

```
1 [metadata]
2 name = plugin-name
3 version = 0.0.1
4 classifiers =
5     Framework :: napari
6
7 [options]
8 packages = find:
9 include_package_data = True
10 install_requires =
11     napari
12     #add all the other packages needed for programming the plugin
13 [options.entry_points]
14 napari.manifest = #tell napari where to find to manifest
15     plugin-name = plugin_name:napari.yaml
```

Once the structure is in place, the plugin can be pip installed and napari will detect and encapsulate the plugin. The plugin will be visible in the plugin's window tab of the napari programme.

## 4 Interactions between plugin and napari

The plugin can only interact with the napari viewer through the *napari.Viewer* class. This class is the napari viewer instance and contains the window, widgets, layers, scroll bars and more. All of which are visible in the viewer. Understanding how to access the different sub-classes and information from this class is the key to coding the plugin interactions. To give the reader a better intuition on how the information is stored, some important sub-classes are made explicit (Code 2). Once the plugin has access, the plugin can also update attributes seen in the viewer. This interaction is bidirectional as plugin functions can connect to viewer *events*. Some of napari's modules, for example *layers*, *Window* and *dims* have events. A module can have several events and each event can be reacted to. Qt implements events with what one

calls signals and slots <sup>3</sup>. The event emits a signal that can then be connected to slot function that executes an operation.

Code 2: Some basic napari.Viewer python access commands

```
napari.Viewer.window #the window widget that contains all dock widgets
napari.Viewer.layers #list of image layers
napari.Viewer.layers[0].source.path #path of the source of first layer
napari.Viewer.layers.events #events class of layers containing all the events a layer can have
napari.Viewer.dims #the napari scroll bar widget
napari.Viewer.dims.events.current_step #the current step event of the scroller
napari.Viewer.dims.events.current_step.connect(function)#connect a slot function to an event
```

## 5 Eda-napari

The plugin programmed for this project is named eda-napari, with eda for event driven acquisition. The whole plugin code can be accessed on GitHub via: <https://github.com/LEB-EPFL/eda-napari>. The main structure of the plugin is based on section 1. The widget code is placed in a separate file `_widget.py`. The structure of the final code is shown in the annex 7.1. The full code of the `_widget.py` file is also available in annex 7.2.

The plugin has two main features: *Frame rate widget* and *Time scroller widget*. Each feature is implemented as a python class that inherits from QWidget. The class is linked to the napari viewer during initialisation by passing as argument `napari_viewer` [8] (Code 3). During the instantiation, napari replaces arguments with values equal to "napari\_viewer" with the viewer class.

The graphical interface is coded in qtpy which is the binding language between Qt (Framework coded in c++) and Python. The main structure consists of creating a Qt layout and then populating it with widgets, for example buttons, labels or figures.

The details of eda-napari's whole class methods will not be explained since the fully annotated code is in the annex and on the Github. However, important implementation choices and references to code extracts will be presented in this section.

The plugin installs automatically napari-aicsimageio [5]. AICSImageIO is an image reader specialised for biological image formats. It is useful for visualising multiple layers in napari. This reader is not necessary for the plugin's functionality. It is also possible to use the builtin reader but it is advantageous to use AICSImageIO. It can be selected as the default reader, when opening a tif file for the first time.

### 5.1 Frame rate widget

Frame rate widget (Code 3) is the major feature of this project. The widget reads metadata from tif files originating from microscope imaging. From a time-lapse image stack, the widget reads the capture time of each individual image. This time data is stored and the frame rate is approximated from this data and then plotted in a dock widget in the napari viewer window. The widget is implemented as a class with an initialisation, class attributes, connections to napari events and class methods.

### Adding a plot to the widget

It is possible to embed a matplotlib canvases in the Qt interface [6]. A Figure canvas is created with axes programmed, as for a typical matplotlib plot. The frame rate is plotted with respect to

<sup>3</sup>Qt Signals and slots documentation: <https://doc.qt.io/qt-6/signalsandslots.html>

time or frame number, with a `QPushButton` allowing to switch between the two. The evolution of the time is also plotted versus the frame number. A vertical line, showing the current frame, is linked to the napari's scroll bar (*dims*) and is updated if the *current\_step* is modified (Code 3, l. 124). For better visibility, a table of the current values is added to the layout under the plots (Figure 3).




Figure 3: Base canvas for the widget plots with interactive buttons and display labels.

## Adding Slow motion shape

A slow motion symbol (Figure 4) is implemented to distinguish two cases: slow frame rate and fast frame rate. Two cases is sufficient since ideally a sample has activity or doesn't. During activity, this is when a higher image frame rate is used. To distinguish these two states otsu thresholding from `skimage`<sup>4</sup> is imported. The package determines the frame rate threshold delimiting the two cases. The easiest way to display an additional image is to use napari viewer's `add_shape()` (Code 3, l. 240) function. This creates a layer automatically and allows quick drawing of polygons. Since several layers can be inserted for one tif file, a `QTimer` (Code 3, l. 92) is necessary to allow the complete upload of all the layers. If the widget is activated before loading an image, the `QTimer` will stop the misplacement of the slow motion layer once data is read.



Figure 4: Slow motion icon.

## 5.2 Time scroller widget

Time scroller widget is the second class widget of this plugin (Figure 5). The widget creates a `QScrollbar` that discretizes time. The goal of the widget is to animate the time-lapse linearly in time. This is useful because napari's scroll bar animates individual frames with a constant time interval without taking into account the actual time from the metadata. The class is coded similarly to the frame rate widget and is set in a second dock widget.

<sup>4</sup>Thresholding documentation: [https://scikit-image.org/docs/stable/auto\\_examples/applications/plot\\_thresholding.html](https://scikit-image.org/docs/stable/auto_examples/applications/plot_thresholding.html)





Figure 5: Base display for the time scroll bar widget.

## Animating linearly in time

For an animation, a relevant time interval needs to be determined. To make this choice, one needs to consider that the animation uses computational power and that the human eye or computer can only detect a limited amount of frames per second. It is therefore unnecessary to overdiscretize. As an arbitrary approximation, 1/4th of the minimum time between image frames is set as the time interval (Code 4, l. 420). A QTimer is linked to a play button and ensures a regular update for the animation. The scroll bar is connected to napari's scroll bar for bidirectional connections. To avoid double calling of events some connections are briefly disconnected (code 5, line: 458-477). For the animation, an index list (Table 1) is created for knowing when to update the viewer's image and not simply just the time incrementation. For speeding up or slowing down the animation, a hybrid method is used. For speeding up the simulation, under a certain critical time interval, *critical\_ms* = 160ms (equivalent to 60Hz computer fps), instead of changing the QTimer's tick time, the step size is increased by a factor 2 (code 5, line: 349-356). The step size is defined as the amount of time intervals skipped after the timer's timeout. Having a tick time under the critical time interval can cause wasteful computations and potential animations problems.

Table 1: Example of frame number indexing in a list for the animation. Between each column, one time interval elapses. The step size is set to one.



## 5.3 Results

The plugin is successfully detected by napari. Both widgets can be opened simultaneously and read data from tif time-lapse images. Example data with variable frame rates were used to verify the execution. The different features of the widgets are demonstrated in the following figures and the AICSImageIO reader is used for it's multiple layer display.

Firstly, figure 6, shows the detection of the plugin in napari. Once the frame rate widget is activated, the plots are displayed on the side, in a dock widget. Since no slow motion icon is visible for the image frame, the threshold frame rate has already been reached. This is the wanted result, since the frame in the viewer is taken during a fast frame rate. The image data used for figure 6 corresponds to a time-lapse of a double channel image sequence. The first channel, in pink, visualises mitochondria fission and the other channel the protein drp1. The camera frame rate used was not implemented with a neural network sensing activity. The frame rate function here is simply for test purposes. On the figure, the two graphs highlight in different ways the rate of image acquisition.






Figure 6: Opening eda-napari Frame rate widget in napari, with mitochondria double channel image sample.

The second widget, the time scroll bar, can also be activated and docked under the main window (Figure 7). The same mitochondria are again visualised but at a different time frame, hence this time a slow motion icon is visible. The x-axes of the frame plot plot has been switched to frame number by clicking on a button.




Figure 7: Frame rate widget and time scroller widget activated with mitochondria time-lapse image and during slow image acquisition.

The final viewer, figure 8, uses test data captured with a periodic evolving frame rate. The regular intervals are clearly seen with the evolution of the side plots. The frame is again in the low frame rate domain, explaining why a slow motion icon appears.

The plugin, as a whole, visualises frame rates from tif files but can also animate linearly in time with the time scroll bar *play >>* button.




Figure 8: Napari viewer with loaded example data sample and periodic varying frame rate.

## 5.4 Limitations and Improvements

The main limitations are due to lengthy computations needed to update the plots. With frame widget activated, napari's scroll bar reactivity is slowed down. It can take up to 1 second for the stop button to stop the animation. The time scroll bar is however not effected. The plugin is overall faster, while using only one single widget.

The animation with the time scroll bar suffers in linear time stepping, when the frame rate widget is used. Between each time tick, many operations need to be executed. During fast animations, the updating of the plots can cause visible halting and disturb the animation.

The plugin is also limited in the format of readable files. More specifically, it only accepts tif files with metadata in s or ms. Additionally, the metadata needs to be structured with specific keywords.

By default, the plugin takes the source of the channel 0 image data and thus assumes that images from other channels have the same acquisition time. This could also be a potential limitation.

General improvements in the plugins reactivity could smoothen the user experience. Depending on the types of microscope data read, additional reading capabilities could also be desired. Currently the plugin has not been published on napari-hub, therefore it can't be downloaded directly on napari but needs to be installed following section 5.7. For publishing in napari, more intensive tests and documentation would be necessary.

## 5.5 Automatic testing

The plugin has been mostly tested with non automatic tests, but ideally intensive automated tests should be implemented. Several automated loading tests have, nevertheless, been implemented with pytest. Pytest operates with fixtures that create an object instance. Functions can then be created to test the instances and simulate user interaction or check for predicted behaviour.

For pytest to evaluate the plugin, it needs access to certain packages. This is why in the *setup.cfg* file of the project there are extra option requirements that are specified for testing. To run the tests, one must simply navigate into the directory *eda-napari* and run in the terminal *pytest*. Pytest detects any test file starting by the name *test\_* in the *eda-napari* folder. For a reader wanting to add tests to the plugin, it can be useful to remember that a *QTimer* postpones the

initialisation of frame rate widget, while data is uploaded. This could induce test failures, if the test does not wait for the widget to load the data.

## 5.6 Issues

A program, especially one with widgets, must be adapted to all sort of different event loops, due to different usages. This is why, for example, the initialisation is programmed with a *try* and *except* statement. This avoids the programme crashing due to unreadable image formats. On the **github** for the project, all current issues are listed and commented. The current detected issues consist of one unwanted interaction between the two sliders and a small loading inconsistency.

## 5.7 Installing the eda napari plugin

This section will be helpful for downloading eda-napari. Eda-napari is not available on napari-hub and must therefore be downloaded manually. The following steps explain the process.

1. First, it is recommended to set up a virtual environment.
2. From **github** download or clone the whole eda-napari folder. All the files must be placed in a folder named “eda-napari” such that the directory structure (Annex 7.1) is respected.
3. Activate the virtual environment.
4. Pip install a backend for Qt. PyQt5 or PySide2 is recommended. The plugin does not specify the backend, but lets the user choose.

```
pip install PyQt5 #or PySide2
```

5. Still within the virtual environment, navigate into the eda-napari directory and install.

```
pip install -e . # "-e" for editing mode
```

6. Napari can now be opened from the terminal and will detect eda-napari.

```
napari
```

## 6 Conclusion

This project illustrates how to create an elaborate npe2 napari plugin, starting from napari’s plugin tutorial. The basics of the Qt frame work and napari’s python classes can be better understood by following the programming of the plugin’s widgets. The plugin is functional and can display, as planned, the frame rate of different time-lapse images. The program can also animate the time-lapse with a separate scroll bar. Theses widgets facilitate the visualisation of event driven acquisition with several means and user interactions. The code could also be used as a template for creating additional npe2 features that handle image metadata.

However, the plugin does still need intensive testing to ensure no bugs can occur and minor computational improvements could be beneficial for the user experience.

## References

- [1] Perkel, J. M. (2021). Python power-up: New image tool visualizes complex data. *Nature*, 600(7888), 347–348. <https://doi.org/10.1038/d41586-021-03628-7>
- [2] Mahecic, D., Stepp, W. L., Zhang, C., Griffié, J., Weigert, M., & Manley, S. (2021). Event-driven acquisition for content-enriched microscopy [Preprint]. *Cell Biology*. <https://doi.org/10.1101/2021.10.04.463102>
- [3] napari contributors (2022). napari: a multi-dimensional image viewer for python. doi:10.5281/zenodo.3555620
- [4] napari contributors (2022). napari tutorial: Your First Plugin. [https://napari.org/plugins/first\\_plugin.html](https://napari.org/plugins/first_plugin.html)
- [5] Allen Institute for Cell Science. (2021). AICSImageIO. <https://github.com/AllenCellModeling/aicsimageio>
- [6] Christopher Nauroth-Kress. (2021), A napari Plugin for visualisation of pixel values over time (t + 3D, t+ 2D) as graphs. [https://github.com/ch-n/napari-time\\_series\\_plotter](https://github.com/ch-n/napari-time_series_plotter)
- [7] Haase, Robert & Fazeli, Elnaz & Legland, David & Doube, Michael & Culley, Siân & Belevich, Ilya & Jokitalo, Eija & Schorb, Martin & Klemm, Anna & Tischer, Christian. (2022). A Hitchhiker’s Guide through the Bio-image Analysis Software Universe.
- [8] napari workshops. (2022). Napari Plugin Accelerator Grant workshop series. <https://chanzuckerberg.github.io/napari-plugin-accel-workshops/>

## 7 Annexes

### 7.1 Project structure




Figure 9: Project structure

### 7.2 The main widget python code

The following code is from `_widget.py` and is accessible also on the [github](#).

#### 1. Imports

```
1 import os
2 import traceback
3 import matplotlib.style as style
4 import matplotlib.pyplot as plt
5 from matplotlib.backends.backend_qt5agg import FigureCanvas #Qt binding not specified (PyQt5,
   PyQt4...)
6 from qtpy.QtWidgets import QWidget, QVBoxLayout, QHBoxLayout, QPushButton, QLabel, QGridLayout
   , QScrollBar
7 from typing import Union
8 import qtpy
9 from qtpy.QtCore import Qt, QTimer
10 import magicgui
11 from magicgui import magic_factory
12 import numpy as np
13 import math
14
15 import napari
16 import tifffile
17 import xmltodict
18
19 from skimage.filters import threshold_otsu
20
21 #Stylesheets
22 #####
23 style.use(str(os.path.dirname(__file__))+'/plot_stylesheet.mplstyle') #get path of parent
   directory of script since plot_stylesheet is in the same directory
24 stylesheet = open(str(os.path.dirname(__file__))+'/q_label_stylesheet.qss', "r")
25 label_style = stylesheet.read()
```

## 2. Frame rate widget

Code 3: Frame rate widget

```
27 #Widgets
28 #####
29
30 #Union is a type: it forms the math union.
31 #It means the widget could be a magicgui widget or a qtpy widget.
32 Widget = Union["magicgui.widgets.Widget", "qtpy.QtWidgets.QWidget"]
33
34 class Frame_rate_Widget(QWidget):
35     """The Frame rate plotter widget.
36
37     This widget is a object inheriting from QWidget. Defining Frame_rate_Widget in the
38     manifest.
39     "napari.yaml" allows Napari to reconise it's existance and display it in the plugin.
40     Upon it's creation the __init__ function ensures
41     the execution of specific functions and plots the frame rates of an OME.tiff file
42     inserted in the napari viewer.
43     """
44
45     def __init__(self, napari_viewer):
46         """Constructor of the Frame_rate_Widget.
47
48         This constructor initialises two blank canvases, the class's viewer to the napari viewer
49         and connects events to functions to allow dynamic plots.
50         A newly inserted file or modification of slider position in napari causes the plot to
51         update. The update after an insertion waits a certain Twait time before executing to
52         ensure
53         napari has time to fully import all the layers. The wait is implemented with a QTimer.
54         """
55         super().__init__()
56         self._viewer = napari_viewer
57         self.image_path=None
58         self.time_data=None
59         self.frame_rate_data=None
60         self.channel=0
61         self.frame_x_axis_time=True
62         self.setStyleSheet(label_style)
63
64         #main plot widget
65         self.layout=QVBoxLayout(self)
66         self.grid_layout = QGridLayout()
67
68         self.button_txt=('Time -> Frame number')
69         self.button_axis_change=QPushButton(self.button_txt)
70         self.button_axis_change.clicked.connect(self.change_axis)
71
72         self.grid_layout.setAlignment(Qt.AlignHCenter)
73         self.grid_layout.setSpacing(2)
74         self.grid_layout.setColumnMinimumWidth(0, 86)
75         self.grid_layout.setColumnStretch(1, 1)
76
77         self.frame_number_label = QLabel("Frame number: ")
78         self.grid_layout.addWidget(self.frame_number_label, 1, 0)
79         self.frame_number_value=QLabel("-")
80         self.grid_layout.addWidget(self.frame_number_value, 1, 1)
81
82         self.frame_time_label = QLabel("Time of frame capture: ")
83         self.grid_layout.addWidget(self.frame_time_label, 2, 0)
84         self.frame_time_value=QLabel("-")
85         self.grid_layout.addWidget(self.frame_time_value, 2, 1)
86
87         self.frame_rate_label = QLabel("Frame rate: ")
88         self.grid_layout.addWidget(self.frame_rate_label, 0, 0)
89         self.frame_rate_value=QLabel("-")
90         self.grid_layout.addWidget(self.frame_rate_value, 0, 1)
91
92         self._init_mpl_widgets()
93         self.layout.addWidget(self.button_axis_change)
94         self.layout.addLayout(self.grid_layout)
95         #QTimer, this Q timer will be used to load data after a certain wait time
96         self.Twait=2500
```

```

92     self.timer=QTimer()
93     self.timer.setInterval(self.Twait)
94     self.timer.setSingleShot(True)#timer only runs once
95     self.timer.timeout.connect(self.init_data)
96
97     #events
98     self._viewer.layers.events.inserted.connect(self.init_after_timer)
99     self._viewer.layers.events.removed.connect(self.update_widget)
100
101 def _init_mpl_widgets(self):
102     """Method to initialise two matplotlib figure canvases with a basic layout and title.
103
104     This method generates a matplotlib FigureCanvas and populates it with a
105     matplotlib.pyplot.figure. The canvas is added to the QWidget Layout afterwards.
106     """
107     self.fig = plt.figure()
108     self.canvas = FigureCanvas(self.fig)
109     self.ax = self.fig.add_subplot(211)
110     self.ax2 = self.fig.add_subplot(212)
111     self.layout.addWidget(self.canvas)
112     self.init_data() #try to init data if it exists
113
114 def init_data(self):
115     try:
116         if self.image_path != self._viewer.layers[0].source.path : #update data if new source
117         is added
118             self.image_path = self._viewer.layers[0].source.path
119             self.time_data=get_times(self)#init times of initial image
120             self.frame_rate_data=self.get_frame_rate()#init frame rate of initial image
121             self.plot_frame_data()
122             self.create_SlowMo_icon()
123             self.slow_mo()
124
125             self._viewer.dims.events.current_step.connect(self.plot_slider_position)
126             self._viewer.dims.events.current_step.connect(self.update_slowMo_icon)
127
128         except (IndexError, AttributeError): # if no image is placed yet then Errors would occur
129         when the source is retrieved
130             print('Meta data not readable')
131         except KeyError:
132             print('Dictionary access in get_times fails. Tif file does not have the adapted keys.')
133     )
134     traceback.print_exc()#prints the info of error
135
136 def init_after_timer(self):
137     self.timer.start(self.Twait) #restarts the timer with a timeout of Twait ms
138
139 def plot_times(self):
140     self.ax.clear() #clear plot before plotting
141     self.ax.set_ylabel('Time [ms]')
142     self.ax.set_xlabel('Frame number')
143     self.ax.set_title('Evolution of elapsed time')
144     self.ax.plot(self.time_data)
145     self.ax.ticklabel_format(axis='y', style='scientific', scilimits=(0,0), useMathText=
146     True')
147     self.line_1=self.ax.axvline(self._viewer.dims.current_step[0],0,1,linewidth=1, color=
148     'indianred')#initilaise a vertical line
149     self.fig.canvas.draw()
150
151 def get_frame_rate(self,unit_frame_r='Hz'):
152     """ Method that returns frame rate.
153
154     Input: unit of frame rate: [kHz] or [Hz]
155     Output: Vector of frame rates in [kHz] or [Hz]
156
157     The frame rate is calculate for each frame. Since the system is discrete, it must be
158     approximated.
159     For the first frame, the frame rate is approximated with the second frame time. For the
160     other frames, the rate is calculated with the previous frame time.""""
161
162     if unit_frame_r=='Hz':
163         conversion_factor =1000 #convert to kHz to Hz
164     elif unit_frame_r=='kHz':
165         conversion_factor=1 #initital data is in kHz since time is displayed in ms
166     else:

```

```

160         print('unit of frame rate not reconised: please use Hz or kHz')
161
162     N_frames= len(self.time_data)
163     frame_rate=[conversion_factor/abs(self.time_data[1]-self.time_data[0])]#the first frame
164     rate
165     for i in range(1,N_frames):
166         frame_rate.append(conversion_factor/abs(self.time_data[i]-self.time_data[i-1]))
167     return frame_rate
168
169 def plot_frame_rate(self,unit_frame_r='Hz'):
170
171     self.ax2.clear()#clear plot before plotting
172     self.ax2.set_ylabel('Frame rate [' +unit_frame_r+']')
173     self.ax2.ticklabel_format(axis='x', style='scientific', scilimits=(0,0), useMathText='
174     True')
175     self.ax2.set_title('Evolution of the Frame rate ')
176     if self.frame_x_axis_time:
177         self.ax2.set_xlabel('Time [ms]')
178         self.ax2.plot(self.time_data,self.frame_rate_data)
179         vline_pos=self.time_data[self._viewer.dims.current_step[0]]
180         txt_text_box='Time of current frame = '+str(self.time_data[self._viewer.dims.
181         current_step[0]])+'[ms]'
182     else:
183         self.ax2.set_xlabel('Frame number')
184         self.ax2.plot(self.frame_rate_data)
185         vline_pos=self._viewer.dims.current_step[0]
186         txt_text_box='Current frame number = '+str(self._viewer.dims.current_step[0])
187
188     self.ax2.set_ylim(0,self.ax2.get_ylim()[1]*1.2)#increase plot for text space
189     self.line_2=self.ax2.axvline(vline_pos,0,1,linewidth=1, color='indianred')#initialise a
190     vertical line
191     props = dict(boxstyle='round', facecolor='wheat', alpha=0.3)
192     self.text_box=self.ax2.text(0.05, 0.95, txt_text_box, transform=self.ax2.transAxes,
193     fontsize=8, verticalalignment='top', bbox=props, color='red')
194     self.fig.canvas.draw()
195
196 def plot_frame_data(self):
197     self.plot_times()
198     self.plot_frame_rate()
199     current_frame=self._viewer.dims.current_step[0]
200     if self.frame_x_axis_time:
201         self.line_2.set_xdata(self.time_data[current_frame])#update according to x axis (time
202         or frame)
203     else:
204         self.line_2.set_xdata(current_frame)
205
206     self.frame_rate_value.setText(str(np.around(self.frame_rate_data[current_frame],5)) + '
207     [Hz]')#init Qlabels
208     self.frame_number_value.setText(str(current_frame))
209     self.frame_time_value.setText(str(self.time_data[current_frame])+' [ms]')
210
211 def plot_slider_position(self,event): #event information stored in "event"
212     """ Method plots and updates slider position on the canvas.
213     Input: event of current_step from slider
214     Output: -
215     After moving the slider on napari viewer, this function is called to update the vertical
216     lines. The vertical
217     lines show the frame rate and capture time of the current image displayed in the napari
218     viewer."""
219     current_frame=event.source.current_step[0]
220     self.line_1.set_xdata(current_frame) #update line
221     if self.frame_x_axis_time:
222         self.line_2.set_xdata(self.time_data[current_frame])#update according to x axis (time
223         or frame)
224     self.text_box.set_text('Time of current frame = '+str(self.time_data[current_frame])+
225     '[ms]')
226     else:
227         self.line_2.set_xdata(current_frame)
228         self.text_box.set_text('Current frame number = '+str(self._viewer.dims.current_step
229         [0]))
230
231     self.frame_rate_value.setText(str(np.around(self.frame_rate_data[current_frame],5)) + '
232     [Hz]')#update Qlabels
233     self.frame_number_value.setText(str(current_frame))

```



```

222     self.frame_time_value.setText(str(self.time_data[current_frame])+ ' [ms] ')
223     self.fig.canvas.draw()
224
225     def change_axis(self):
226         self.frame_x_axis_time=not self.frame_x_axis_time
227         if self.frame_x_axis_time:
228             button_txt='Time -> Frame number'
229             self.current_frame_txt='Current frame number = '
230
231         else:
232             button_txt='Frame number -> Time'
233         self.button_axis_change.setText(button_txt)
234         self.plot_frame_rate()
235
236     def create_SlowMo_icon(self):
237         triangle=np.array([[20, 60], [60, 60], [40, 90]])
238         rectangle=np.array([[20, 40],[60, 40],[60,25],[20,25]])
239         polygon=[triangle,rectangle]
240         self._viewer.add_shapes(polygon, shape_type='polygon', face_color='white',edge_width=2,
241                                edge_color='black', name='Slow motion')
242         self.slow_mo_channel=self._viewer.layers.index('Slow motion')
243         self._viewer.layers[self.slow_mo_channel].visible=False #init to invisible
244
245     def update_slowMo_icon(self,event):
246         if(self.slow_mo_array[event.source.current_step[0]]):
247             self._viewer.layers[self.slow_mo_channel].visible=True
248         else:
249             self._viewer.layers[self.slow_mo_channel].visible=False
250
251     def slow_mo(self):
252         size=len(self.time_data)
253         self.slow_mo_array=np.empty(size)
254         thresh=threshold_otsu(np.array(self.frame_rate_data))#threshold to determine weather
255         slow motion or fast speed
256         for i in range(0,size):
257             if self.frame_rate_data[i]>=thresh:
258                 self.slow_mo_array[i]=False
259             else:
260                 self.slow_mo_array[i]=True
261
262     def update_widget(self,event):
263         """ Method updates the widget in case of layer deletion.
264         Input: event created my deleted layer.
265         Output: -
266         If slow motion is removed then the napari viewer is disconnected to the slow motion
267         shape.
268         If all layers are removed the dock plugin is removed.
269         """
270         if not ('Slow motion ' in event.source):
271             self._viewer.dims.events.current_step.disconnect(self.update_slowMo_icon)
272         if len(event.source)==0:
273             try:
274                 self._viewer.window.remove_dock_widget(self)
275                 self.image_path=None
276             except:
277                 print('Dock already deleted')

```

### 3. Time scroller widget

Code 4: Time scroller widget

```

278 class Time_scroller_widget(QWidget):
279     """Time_scroller_widget class is a widget that creates a time scroll bar. The scroll bar
280     allows to animate stacks of
281     images linearly with time. Similary to the napari scroll bar, it has a play, stop, next and
282     previous button."""
283
284     def __init__(self, napari_viewer):
285         """Constructor of the Time_scroller_widget.
286         This constructor initialises the button widgets and scroll bar widget in a QHBoxLayout.
287         It also initialise the time of the frames
288         from the image data available. Some signals are also defined to allow interaction and
289         automatic updates between the current layer and

```

```

287     the Time_scroller_widget widget. A QTimer is defined to controll the animation.
288     """
289     super().__init__()
290     self._viewer = napari_viewer
291     self.image_path=None #folder path of the image data
292     self.time_data=None
293     self.channel=0
294     self.number_frames=None
295     self.show_time=50#animation default display time ms
296     self.time_interval=None # time of the discretised time interval [ms]
297     self.interval_frames_index=[]#discretised time interval filled with an index the the
    different frames.
298     self.step=1
299     self.critical_ms=160
300     self.timer=QTimer(self)
301     self.timer.timeout.connect(self.play_step)
302
303     self.layout=QHBoxLayout(self)
304     self.setMinimumWidth(500)
305     self.setMaximumHeight(100)
306
307     self.create_bottom_dock_button()
308     self.create_slow_down()
309     self.create_play_button()
310     self.create_speed_up()
311     self.create_time_scroller()
312     self.create_axis_label()
313
314     self.layout.addWidget(self.bottom_dock_button)
315     self.layout.addWidget(self.slow_down_button)
316     self.layout.addWidget(self.play_button)
317     self.layout.addWidget(self.speed_up_button)
318     self.layout.addWidget(self.time_scroller)
319     self.layout.addWidget(self.axis_label1)
320     self.layout.addWidget(self.axis_label2)
321
322     self.init_data()
323     self._viewer.layers.events.inserted.connect(self.init_data) #init data when layer is
    inserted
324
325     def create_bottom_dock_button(self):
326         self.bottom_dock_button=QPushButton('')
327         self.bottom_dock_button.setMaximumWidth(35)
328
329     def move_dock_to_bottom(self):
330         self.parentWidget().parentWidget().addDockWidget(Qt.BottomDockWidgetArea, self.
    parentWidget()) # init position in QDockWidget (parent) to bottom in QWindow
331         self.bottom_dock_button.deleteLater()
332
333     def create_slow_down(self):
334         self.slow_down_button=QPushButton('x 0.5')
335         self.slow_down_button.setMaximumWidth(50)
336
337     def slow_animation(self):
338         if self.step == 1:
339             self.show_time = self.show_time*2
340             if self.play_button_txt=='Stop': #if program is playing
341                 self.timer.stop()
342                 self.timer.start(self.show_time) #set new show_time
343         else:
344             self.step=self.step/2
345     def create_speed_up(self):
346         self.speed_up_button=QPushButton('x 2')
347         self.speed_up_button.setMaximumWidth(50)
348
349     def speed_animation(self):
350         if self.show_time >=self.critical_ms:
351             self.show_time = self.show_time*0.5
352             if self.play_button_txt=='Stop':
353                 self.timer.stop()
354                 self.timer.start(self.show_time)
355         else:
356             self.step=self.step*2
357
358     def create_time_scroller(self):

```

```

359     self.time_scroller= QScrollBar(Qt.Horizontal)
360     self.time_scroller.setMinimum(0)
361     self.time_scroller.setSingleStep(1)
362     self.time_scroller.setMinimumWidth(150)
363
364     def create_play_button(self):
365         self.play_button_txt='Play >>'
366         self.play_button=QPushButton(self.play_button_txt)
367         self.play_button.setMaximumWidth(60)
368
369     def create_axis_label(self):
370         self.axis_label1=QLabel('End time')
371         self.axis_label2=QLabel('Current time')
372
373     def init_data(self):
374         """This method initialises all the additional data after an image stack is available and
375         readable.
376         """
377         try:
378             if self.image_path!=self._viewer.layers[0].source.path: #Only inits data if the layer
379             is new
380                 self.image_path = self._viewer.layers[0].source.path
381                 self.time_data=get_times(self)#init times of image stack
382                 self.number_frames=len(self.time_data)
383                 self.init_time_interval()
384                 self.set_frames_index()
385                 #time scroller
386                 self.time_scroller.setMaximum(len(self.interval_frames_index)-1)
387                 idx=self.interval_frames_index.index(self._viewer.dims.current_step[0])
388                 self.time_scroller.setValue(idx) #set init position of scroller
389                 #init label2
390                 self.axis_label2.setText('| '+str(self.time_data[-1])+ ' [ms] ')
391                 width2 = self.axis_label2.fontMetrics().boundingRect(self.axis_label2.text()).
392                 width() #max width of text
393                 self.axis_label2.setFixedWidth(int(1.1*width2))
394                 #init Label1
395                 self.axis_label1.setText(str(self.time_scroller.value()*self.time_interval))
396                 self.axis_label1.setFixedWidth(int(0.7*width2))
397
398                 #events
399                 self.bottom_dock_button.clicked.connect(self.move_dock_to_bottom)
400                 self.slow_down_button.clicked.connect(self.slow_animation)
401                 self.play_button.clicked.connect(self.play)
402                 self.speed_up_button.clicked.connect(self.speed_animation)
403                 self._viewer.dims.events.current_step.connect(self.update_scroller_from_dims)#
404                 link window srolle to time srolle
405                 self.time_scroller.valueChanged.connect(self.update_scroller_from_scroller) #
406                 link
407                 self._viewer.layers.events.removed.connect(self.update_widget)
408                 self.data_is_available=True
409
410         except (IndexError, AttributeError): # if no image is found then an index Error would
411         occur
412             print('Meta data not readable')
413         except KeyError:
414             print('Dictionary access in get_times fails. Tif file does not have the adapted keys.
415             ')
416         traceback.print_exc()#prints the info of error
417
418     def init_time_interval(self):
419         """This method sets the time discretisation interval of the system, for the animation
420         and scroll bar.
421         self.time_interval is initiliased as the 1/4 of the minimum time between to images
422         frames. This creates a revelant discretisation of time with
423         respect to the frame rates.
424         """
425         diff=[]
426         for i in range(1,self.number_frames):
427             diff.append(self.time_data[i]-self.time_data[i-1])
428         min_diff=min(diff)
429         self.time_interval = math.floor(min_diff/4)#this makes sure a relevant discretisation of
430         time is made for the animation
431
432     def set_frames_index(self):
433         """This method sets self.interval_frames_index with image frames numbers. Each index

```

```

424     corresponds the appropriate image frame
425     that should be displayed in the discretized time.
426     """
427     frame_index=0
428     t=0
429     self.interval_frames_index=[0]
430     while frame_index < self.number_frames-1:
431         t+=self.time_interval
432         if self.time_data[frame_index+1] < t:
433             frame_index+=1
434         self.interval_frames_index.append(frame_index)
435
436 def play_step(self):
437     """This method advances the time of 1 time interval and takes care of updating the
438     current displayed frame if necessary.
439     """
440     if self._viewer.dims.current_step[0]== self.number_frames-1:
441         self._viewer.dims.set_current_step(0, 0)#restart at frame 0
442         self.time_scroller.setValue(0)
443     else:
444         self.time_scroller.setValue(self.time_scroller.value()+self.step)
445         if self.interval_frames_index[self.time_scroller.value()] != self._viewer.dims.
446         current_step[0]: #update viewer
447             self._viewer.dims.set_current_step(0, self.interval_frames_index[self.
448             time_scroller.value()])
449
450 def play(self): # play and stop method
451     if self.play_button_txt == 'Play >>':
452         self.timer.start(self.show_time)
453         self.play_button_txt = 'Stop'
454         self.play_button.setText(self.play_button_txt)
455         self.play_button.setMaximumWidth(80)
456
457     else: #stop
458         self.timer.stop()
459         self.play_button_txt = 'Play >>'
460         self.play_button.setText(self.play_button_txt)
461
462 def update_scroller_from_dims(self):
463     idx=self.interval_frames_index.index(self._viewer.dims.current_step[0])
464     self.time_scroller.valueChanged.disconnect(self.update_scroller_from_scroller)#avoid
465     double calling
466     self.time_scroller.setValue(idx)
467     if self._viewer.dims.current_step[0]== self.number_frames-1:
468         self.axis_label1.setText(str(self.time_data[-1]))
469     else:
470         self.axis_label1.setText(str(self.time_scroller.value()*self.time_interval))
471     self.time_scroller.valueChanged.connect(self.update_scroller_from_scroller) #reconnect
472
473 def update_scroller_from_scroller(self):
474     if self.interval_frames_index[self.time_scroller.value()] != self._viewer.dims.
475     current_step[0]: #update viewer #one
476         self._viewer.dims.events.current_step.disconnect(self.update_scroller_from_dims) #
477         avoid double calling
478         self._viewer.dims.set_current_step(0, self.interval_frames_index[self.time_scroller.
479         value()])
480         self._viewer.dims.events.current_step.connect(self.update_scroller_from_dims) #
481         reconnect
482
483     if self._viewer.dims.current_step[0]== self.number_frames-1:
484         self.axis_label1.setText(str(self.time_data[-1]))
485     else:
486         self.axis_label1.setText(str(self.time_scroller.value()*self.time_interval))
487
488 def update_widget(self, event):
489     if len(event.source)==0:
490         try:
491             self._viewer.window.remove_dock_widget(self)
492             self.image_path=None
493         except:
494             print('Dock already deleted')

```

#### 4. Read tiff file metadata

To read the meta data from tiff files, the following function has been implemented. The meta

data is read with TiffFile and the xltdict package. With the correct dictionary keys the time data of the file is then accessed. This function is not a class method and is accessible from both widgets.

Code 5: Get\_times metadata access function

```

488 def get_times(widget):
489     """Method that gets the capture time from the metadata.
490
491     Input:—
492     Output: Vector of time metadata [ms] of images found at given image path and channel.
493
494     The times of each image stack from a ome.tif file is read in [ms] or [s] and then returned
495     in [ms].
496     The times are taken from a given channel. The data can only be read from an ome.tif file.
497     The Offset
498     from time t=0 subtracted to the times before it is returned.
499     The following code is inspired from the solution for reading tiff files metadata from Willi
500     Stepp.
501     """
502     times=[]
503     with tiffFile.TiffFile(widget.image_path) as tif:
504         XML_metadata= tif.ome_metadata #returns a reference to a function that accesses the
505         metadata as a OME XML file
506         dict_metadata=xltdict.parse(XML_metadata) #converts the xml to a dictionary to be
507         readable
508         num_pages=len(tif.pages) #the number of images stacked
509         for frame in range(0,num_pages):
510             #time should be in either s or ms
511             if float(dict_metadata['OME']['Image']['Pixels']['Plane'][frame]['@TheC'])==widget.
512             channel: #checks if correct channel
513                 frame_time_unit=dict_metadata['OME']['Image']['Pixels']['Plane'][frame]['
514                 @DeltaTUnit']
515                 if frame_time_unit== 's' :
516                     convert_unit_to_ms=1000
517                     times.append(convert_unit_to_ms*float(dict_metadata['OME']['Image']['Pixels']['
518                     Plane'][frame]['@DeltaT']))
519                 elif frame_time_unit == 'ms':
520                     convert_unit_to_ms=1
521                     times.append(convert_unit_to_ms*float(dict_metadata['OME']['Image']['Pixels']['
522                     Plane'][frame]['@DeltaT']))
523                 else:
524                     print('Time units not in ms or s but in '+ frame_time_unit+'. A conversion to
525                     ms or s must be done.')
526
527     times = [x - times[0] for x in times] #remove any offset from time=0
528     return times

```

## 7.3 Debugging

The plugin was coded with the python editor Virtual Studio Code. VS code allows connecting with git and debugging. Setting up a *debug.py* file was useful for gaining time during tests. The debug file opens up napari with the debugger and loads automatically an image. With the debugger the plugin runs much slower than directly being opened in the terminal. All the variables and plugin objects can be accessed any time during the code execution, which helps for locating errors.